



Trampoline Trimilin Vario Plus 120cm-29cm



With **VARIO system** for individual tension of the jumping mat. With these rubber cable models, you can change the tension of the rubber cable trampolines yourself in just a few simple steps. Mini trampolines with the patented Vario system are individually adjustable according to body weight. Tighter or softer, for different applications. With the Trimilin Vario 120 - **with retractable feet** - the basic setting of the suspension (degree of softness) can be changed in both directions, i.e. tighter or softer. This allows you to adjust your Trimilin-Vario individually, e.g. if several people with very different body weights want to train on the trampoline. All without replacing cables or rubber rings! The Trimilin-Vario trampoline is ideal for dynamic swinging and jumping and suitable for the whole family. From the very young to adults. Combines soft suspension with active-dynamic rebound and high fun factor. Advantageous for exercises while lying down and with a partner is the large diameter of 120cm.

CHF 479.00

VARIO system - VARIO is a new type of suspension and tensioning system with 12 individual rubber cables. This allows you to individually adjust the tension of the jumping mat tighter or softer, without having to change cables or rubber rings.

Training on the mini trampoline

- helps to boost the fat burning process
- enables a gentle start to losing weight
- helps to increase the personal fitness and condition individually
- for noticeable and visible success with gentle training



- offers motivation through the fun of swinging
- improves coordination and balance
- in the space-saving "mini home gym"

Healthy fat burning through gentle exercise

Exercise is an important factor in the treatment of overweight and obesity. Obesity is defined as an increase in body fat above normal levels. According to a 2012 study by the Robert Koch Institute, about 67% of men and 53% of women in Germany are overweight or obese. Obesity has increased worldwide over the past two decades, especially among men and young adults. Increasing physical activity is an important item in the energy balance. In particular, endurance exercise - done consistently over months and years - serves to reduce weight.

Trampoline with particularly joint-gentle suspension

Gentle swinging on the Trimilin trampoline is also used in osteoporosis prevention and to alleviate arthritis. The special suspension is gentle but not too soft and therefore provides the necessary safety and stability for all therapeutic applications. This mini trampoline is ideally suited for use in physiotherapy, rehabilitation and physical therapy. It serves to stabilize ankle, knee and hip joints and supports the targeted development of atrophied muscles after operations.

And by the way: anti-aging on the mini-trampoline

A welcome side effect of training on the mini-trampoline is fat reduction and body shaping. The stimulation of cells, lymph and metabolism while swinging and jumping (rebounding) on the mini trampoline has a rejuvenating effect on skin and appearance. The tissue tightens. The general posture improves due to the strengthening of the basic and supporting muscles. Especially on the legs this effect is quickly noticeable, even cellulite can be reduced.

Balance and coordination through regular trampoline swinging

Trampoline training on the Trimilin-med also trains the sense of balance and gently increases mobility. Every cell is stimulated by the ups and downs on the elastic jumping mat, the entire metabolism, lymph and heart circulation are stimulated. Swinging on the mini trampoline sets the whole body in harmonic vibration and helps the organism to dissolve blockages.

Increase cardio fitness with the trampoline

Gentle cardio and fitness training on the trampoline helps to support the venous pump and relieve the heart. Varicose veins also respond successfully to the operation of the venous pump. Especially for older people, trampoline swinging provides relief from joint stiffness, digestive problems, vascular and cardiac problems, too low or too high blood pressure.

Exercise on the trampoline - simple but effective

Gentle or athletic - exercise on the mini trampoline is possible at any age and is very effective. The increased muscle tone improves one's body awareness and physical and mental agility. A few minutes of Trimilin training every day will have an effect on fitness and charisma in a short time.

- Endurance and fitness training
- Physiotherapy and rehabilitation
- Especially for light persons
- Walking and running training possible
- Gentle on intervertebral discs and joints
- Gentle rebounding



Short info rubber cable:

- soft
- longer spring travel- slower rhythm
- for special sensitivity of spine and joints
- also suitable for very low body weight

Advantages of rubber cables:

Rubber is softer than steel. Therefore it has a longer spring travel and a slower oscillation frequency. This is particularly favorable for slow, finely dosed vibration. For example, if joints or spine are very sensitive to pain. The lower the body weight, the more recommended this type of suspension is. Especially for small children or light adults.

For Trimilin are used specially sheathed, particularly durable rubber ropes. High stretchability, elasticity and long durability are the outstanding quality features of Trimilin rubber cables. Conventional rubber rings are surrounded only with a thin textile layer. Each Trimilin rubber cable, on the other hand, is encased in a highly elastic synthetic fabric that extends its service life many times over. This ensures the outstanding jumping characteristics and the active-dynamic rebound typical of Trimilin with a particularly high fun factor.

- with elastic rubber cable
- with VARIO system - tension of the jumping mat can be adjusted tighter or softer depending on the application and body weight
- **with foldable feet** for space-saving storage
- diameter 120cm
- diameter of jumping mat 85cm
- recommended up to a body weight of approx. 120kg
- leg length 20cm
- small distance to the ground of approx. 29cm
- available with grab handle
- degree of softness 1-6
- weight 10kg
- jumping mat: black or blue
- border cover: silver - other colors for the outer border available for an extra charge
- incl. tensioning instructions and tensioning tool