



Tanita DC 360 P body fat scale



With the help of the clinically accurate Tanita Dual Frequency BIA technology, a complete body analysis can be carried out in 15 seconds.

CHF 3'750.00

The DC 360 offers total body measurements in just a few seconds. The easy-to-use console is designed to guide the user safely through the individual steps of data entry. This makes it ideal for frequent use. The printout allows for quick analysis of key results so that the user can understand their health status at a glance. Measurement results such as basal metabolic rate and metabolic age have proven to be very motivating. All measurements can be printed using the integrated thermal printer or automatically transferred to software for data logging and long-term analysis.

Total body measurements:

- Body fat
- Fat mass
- Fat-free mass
- Muscle mass
- Total body water %
- Body Mass Index
- Bone mass
- Body composition value
- Display visceral fat
- Basal metabolic rate
- Basal metabolic rate indicator
- Metabolic age

Product features:

- Whole body composition analysis in 15 seconds



with clinically accurate Tanita dual frequency BIA technology

- Built-in printer provides instant readout of results with header analysis.
- Results are automatically saved to the SD card , sent to a PC or printed.
- GMon Pro software compatible, enabling trend analysis, health risk assessment and full data management.
- Low profile platform for added stability
- Max. Weight capacity 270kg with 100g precision
- 5 year warranty